












ACCESSORIES AND PARTS

	Scale	HMP-S200 model. Calibrated from 0~600 lbs and will fit all HMP standard floor lifts. Accurate digital read outs in either pounds or kilograms. Assembly hardware included with complete instruction.
	Batteries	HMP-4005ET-BAT Replacement Batteries (2 required, fits all HMP lifts)
	Charger	HMP-4005ET - ADR Economic Charger (fits all HMP lifts)
	Hand Control	HMP-4005ET - HC Economic Hand Control (fits all HMP lifts)
	Control Box	HMP-4005ET - BOX 400 lbs. Lift Control Box with PCB
	Control Box for SL600	HMP-SL600 - BOX 600 lbs. Stand Assist Lift Control Box with PCB
	Actuator	LA31-PLT For 400 lbs. Patient Lifts, item # :HMP-4004H ; HMP-4004ET ; HMP-4005E
	Actuator	LA31-PLD For 300 lbs. Patient Lifts, item # :HMP-4004H ; HMP-4004ED
	Actuator	LA31-SL For 400 lbs. Stand Assist Lift, item # :HMP-SL400
	Actuator	LA34-PL For 600 lbs. Patient Lift, item # :HMP-4006E
	Actuator	LA34-SL For 600 lbs. Stand Assist Lift, item # :HMP-SL600

	HMP-4004H	HMP-4004ET	HMP-4004ED	HMP-4005E	HMP-4006E BARIATRIC	HMP-SL400	HMP-SL600 BARIATRIC
Scale	✓	✓	✓	✓	✓		
Batteries	✓	✓		✓	✓	✓	✓
Charger	✓	✓		✓	✓	✓	✓
Hand Control	✓	✓		✓	✓	✓	✓
Control Box	✓	✓		✓	✓	✓	
Control Box for SL600							✓
Actuator LA31-PLT	✓	✓		✓			
Actuator LA31-PLD	✓		✓				
Actuator LA31-SL						✓	
Actuator LA34-PL					✓		
Actuator LA34-SL							✓