



Home Medical Products Inc.

# HANDI BELT



## HMP-TS307

### HMP Handi Belt

is a padded belt with a non slip inner surface. Machine washable , the comfortable to grip, uniquely positioned handles provide a secure and easy way to perform multiple manual transfers.

HMP Handi Belt has uniquely placed handles which are designed to minimize caregiver wrist sprain by keeping the wrists in a straight position. Available in two sizes.



Item No.	SIZE	MEASUREMENTS
Ts30720	Standard	5" x40" adjustable 29" -45"
Ts30730	Large	5" x40" adjustable 35" -55"



Home Medical Products Inc.

# TRANSFER SHEET/QUILT



**HMP-TS301**

A continuous loop of low friction fabric designed to significantly reduce the force required to reposition and transfer patients.

Specially woven fabric is made to glide easily in all directions.


Multi direction **HMP Slide Sheets** are available in number of convenient sizes. Typically the slides is provided as a basic nylon treated sheet. For added comfort HMP Slide Sheets are also available in padded quilt versions.

HMP Slide Sheets are easy to use and clean and provide an unsophisticated, inexpensive but much needed aid to improving daily nursing and care chores.



**HMP-TS302**

Three different sizes



Item	Type	SIZE	MEASUREMENTS
Ts30110	Sheet	Small	28.5" x 27.5"
Ts30120	Sheet	Medium	42.5" x 31.5"
Ts30130	Sheet	Large	53.0" x 39.5"
Ts30210	Quilt	Small	23.6" x 23.6"
Ts30220	Quilt	Medium	23.6" x 55.1"
Ts30230	Quilt	Large	23.6" x 74.8"

# SEAT/STAND TURN



## HMP-TS30810

### SEAT TURN


STANDARD Top is quilted cotton.  
Bottom is non-slip material.  
Inside is low friction material.  
One size: 15.75" (40cm) diameter



## HMP-TS30820

### STAND TURN

STANDARD Top and bottom surfaces  
are non-slip material.  
One size: 15.75" (40cm) diameter



## HMP-TS30830

### SWIVEL CUSHION

STANDARD Top and bottom with sponge  
covered by quilt cotton.  
Inside is swivel balls easy to transfer.  
One size: 15.75" (40cm) diameter

