



Home Medical Products Inc.

SHOWER CHAIRS



HMP-5001

W/O Backrest



HMP-5002

W/ Backrest



HMP-5003

W/ Small Backrest



HMP-5004

U Shape w/o Backrest



HMP-5005

U Shape w/ Backrest

"Quick Release" design available.



Home Medical Products Inc.

SHOWER BENCH



HMP-5007

Maxi-Drain Transfer Bench

Weight Bearing Capacity: 250lbs

Flared legs enhance stability.

Surface Processing : Satin finish

Seat Size:(16" X 8.5") X 3; Backrest Size: 16.5" X 7"



HMP-5009

Imperial Transfer Bench

Weight Bearing Capacity: 420 lbs

Flared legs and tips of suckers enhance stability.

Surface Processing : Satin finish

Seat Size: 30" X 16.5"
Backrest Size: 16.5" X 7"



HMP-5010

Imperial Transfer Bench

Weight Bearing Capacity: 250lbs

Flared legs and tips of suckers enhance stability.

Surface Processing: Satin finish

Seat Size: (16" X 7-3/4") X 3 ;
Backrest Size: 16" X 7-3/4"



HMP-5012

U Shape Transfer Bench

Bearing Capacity: 300 lbs.

Flared legs and tips of suckers enhance stability.

Surface Processing : Satin finish

Seat Size: (16" X 7-3/4") X 1 and U-Shape seat;
Backrest Size: 16" X 7-3/4"



Home Medical Products Inc.

SHOWER CHAIR



HMP-5006

Deluxe Padded Bath Chair

Weight Bearing Capacity: 250 lbs

Flared legs enhance stability.

Surface Processing :
Satin finish

Adjustable height : 1"
increments from 17" -
21" (ground to seat)



HMP-5009

Deluxe Padded Bath & Commode Chair

Weight Bearing Capacity: 250lbs

Flared legs enhance stability.

Surface Processing :
Satin finish

Seat Size: (16" X 8.5")
X 3 ; Backrest Size :
16.5" X 7"



HMP-5011

Round Seat Shower Chair

Weight Bearing Capacity: 250 lbs

Flared legs enhance stability.

Surface Processing :
Shiny or Satin finish
available

Seat Size : Diameter
12.5"

Adjustable height: 1"
increments from 16"-
22" (ground to seat)

"Tool Free" assembly available.

