

# STANDER:

- Positioning Chairs.....47,48
- Bolster Chair.....49
- Tilt Positioning Desk.....50
- Vertical/Prone/Supine  
Standers .....51~54



## Positioning Chair

## GIRAFFE

## HMP-WP001

## STANDER



Push Bar



- Standard Supply :**
- Tray and Hand Bar
  - Safety Belt
  - Lateral support
  - Leg separator

- Optional Items :**
- Heel Loop of Chest
  - Heel Loop of Waist

### FEATURES:

- Material : Finnish forest Birch Plywood , 17mm , 13ply
- The Chair can be Adjusted forward or backward freely
- The Depth of Seat is Adjustable
- The Width of Seat is Adjustable from Lateral Supports
- The Height of Seat is Height Adjustable
- The Length & Angel of Footplate is Adjustable
- The Reclining Angle of Backrest Has Four Selections. (-5 , 0 , 5 , 12 degree)
- The Angel of Seat is Adjustable by Four Selections. (-5 , 0 , 5 , 12 degree)
- The Height of Armrest is Adjustable
- The Height of Seat Back is Adjustable
- You can move the Tray & Hand Bar Forward or Backward
- Custom-Made Special Size On Request
- Comply with BS-4875 Part 1:2001 Requirement



Safety Belt

Item	Size 1	Size 2	Size 3
Seat Depth	8" -14"	12" -16"	14" -18"
Seat Back Height	10" -16"	14" -20"	16" -20"
Seat Width	8" -14"	12" -16"	14" -18"
Seat Height	8" -14"	12" -18"	16" -20"



Footplate & Heel Loop



Tray and Hand Bar



Head Support



Lateral Support



Seat Pad



Leg Separator

## STANDER

### FEATURES:

- Material : Finnish forest Birch Plywood , 17mm , 13ply
- The Chair can be adjusted forward or backward freely
- The Depth of Seat is Adjustable
- The Width of Seat is Adjustable from Lateral Supports
- The Height of Seat is Height Adjustable
- The Height & Angel of Footplate is Adjustable
- The Reclining Angle of Backrest is Adjustable from 90 degree to 135 degree
- The Angle of Tilt-in is Adjustable
- The Height of Armrest is Adjustable
- The Height of Back is Adjustable
- You can move the Tray & Hand Bar Forward or Backward
- Custom-Made Special Size On Request
- Comply with BS-4875 Part 1:2001 Requirement

## HMP-WP002



## HIPPO

- Standard Supply :
- Tray and Hand Bar
  - Safety Belt
  - Lateral support
  - Leg separator
  - Head Support
- Optional Items :
- Heel Loop of Chest
  - Heel Loop of Waist

## Positioning Chair



Reclining & Tilt-in

Item	Size 1	Size 2	Size 3
Seat Depth	8" -14"	12" -16"	14" -18"
Seat Back Height	10" -16"	14" -20"	16" -20"
Seat Width	8" -14"	12" -16"	14" -18"
Seat Height	8" -14"	12" -18"	16" -20"



Safety Belt & Lateral Support



Footplate & Heel Loop



Tray and Hand Bar



Depth of Seat is Adjustable



Head Support



Seat Pad



Leg Separator

## Bolster Chair

## RABBIT

## HMP-WP003

## STANDER



The Seat Height and Back Height are Adjustable



### FEATURES:

- Material : Finnish forest Birch Plywood , 17mm , 13ply
- Tray Dimension : 65 cm × 82 cm
- Height Adjustable Backrest
- Height Adjustable Seat
- Height Adjustable Tray
- Adjustable footrest
- Moving Wheels
- Custom-Made Special Size On Request
- Comply with BS-4875 Part 1: 2001 Requirement



Tray

Item	Size 1	Size 2	Size 3
Roller Length	60cm	60 cm	60 cm
Roller Diameter	20 cm	25 cm	30 cm
Roller Height	8" -14"	10" -16"	12" -18"
Hips Width	6" -14"	10" -16"	12" -18"



Moving Wheels



Footplate



Hip Support Pads



Roller

## STANDER

HMP-WP004

Lily

Tilt Positioning Desk

### FEATURES:

- Material : Finnish forest Birch Plywood , 17mm , 13ply
- Range of Height: 47~77cm (w/o wheels)
- Wheel Height: 15 cm
- Dimension of Table: 18×60 cm
- Inner width of Table: 65 cm
- Inner depth of Table: 50 cm
- Table tilt Degree: 0~40 degree
- Inner arc Degree: 45×20 cm
- The Height and Depth of the Book-Box is Adjustable
- Custom-Made Special Size On Request
- Comply with BS-4875 Part 1:2001 Requirement



The Angel of Tray is Adjustable



Moving Wheels



The Height of Tray is Adjustable



Book-Box

## Vertical Stander

## AQUARIUS

## STANDER

### FEATURES:

- Material : Finnish forest Birch Plywood , 17mm , 13ply
- Wide Height Range of Hip and Chest Support
- Quick Release Lock
- Easy Moving Wheels
- Hip Width Range ( 8" -17" )
- Total Length 84cm
- Tray : 51cm×55cm ( Height & Depth Adjustable )
- Custom-Made Special Size On Request
- Comply with BS-4875 Part 1:2001 Requirement

### Option :

- Chest Supports
- Tray



Item	Size 1	Size 2
Heel Support Height	5" -13"	9" -21"
Hip Support Height	9" -23"	13" -31"
Chest Support Height	13" -31"	21" -41"
Tray Height	60-80 cm	80-120 cm
Max. User Height	80-120 cm	110-160 cm
Base Total Width	55 cm	60 cm



Tray



Heel Stop



Moving Wheels



Quick Release Lock



Hip Support Pads and Chest Support Pads



Knee Support Pads

## STANDER

## HMP-WP005

## KANGAROO

## Prone Stander

### FEATURES:

- Material : Finnish forest Birch Plywood , 17mm , 13ply
- Standing Angel is Adjustable from 90 degree to 120 degree
- Independent Knee Support Pad and its Width is Adjustable
- The Width & Height of Hip Support Pad and Chest Support Pad are Adjustable
- Shoe Shape of Footplates. The Distance between two footplates can be Adjusted
- The Angel of Footplate is Adjustable
- Tray Dimension : 60cm×60cm
- The Angel of Tray is Adjustable
- Kneel Posture Plate (Option)
- Custom-Made Special Size On Request
- Comply with BS-4875 Part 1:2001 Requirement



The Angel of Tray is Adjustable

Item	Size 1	Size 2	Size 3
Tray Height	60-80 cm	80-105 cm	100-125 cm
Base Length	90 cm	100 cm	110 cm
Base Width	60 cm	65 cm	70 cm
Height for Children	90-120 cm	110-140 cm	125-170 cm



Kneel Posture Plate



Moving Wheels



Footplate



Angle Adjustable Support



Chest Support Pads




Hip Support Pads



Knee Support Pads

## Supine Stander



Angle Adjustable Support



Chest Pads & Hip Pads

## EAGLE



## STANDER

### FEATURES:

- Material : Finnish forest Birch Plywood, 17mm , 13ply
- Standard Back Board Height : 120cm
- Back Board Width : 56cm
- Extension Board : 40cm
- Tilt Angle: 3~35 Degree
- Heel Height Range : 5" -20"
- Trunk Support Width : 6.5" ~18"
- Angle Adjustable Tray
- Tray Dimension : 70cm×55cm
- Custom-Made Special Size On Request
- Comply with BS-4875
- Part 1:2001 Requirement

### Standard Supply :

- Knee Pads or knee Support Pads
  - Chest Pads, Hip Pads
  - Angle Adjustable Board, Footrest
- Optional Items :
- Head Extension Board 40cm )
  - Head Pads
  - Tray



Head Extension Board And Head Pads



Knee Support Pads



Knee Pads



Angle Adjustable Tray



Angle Adjustable Board



Footrest



## STANDER

## HMP-SP010

## DIY Colorful Vertical Stander

### FEATURES:

- Colorful painting for visual attraction
- DIY design
- Soft Chest Belt & Hip Belt
- Easy to adjust forward, backward, upper and lower position of Quick Release Knee Stopper
- 4 wheels with brake
- Option : 360° Heel stopper with straps.

### Standard Supply :

- Detachable Armrest
  - Tray
  - Chest Belt and Hip Belt
  - Quick Release Knee Stopper
- 
- Option Item : Quick Release Knee Stopper



Size	Frame Width	Tray Height	User Height
1	40cm	60-90cm	90-130cm
2	50cm	90-130cm	130-170cm



Tray



Quick Release Knee Stopper



Detachable Armrest



Chest Belt and Hip Belt



360° Heel Stopper and Straps