

# Suction Bathtub Handle



**"RED"** is a warning signal request re-mount again



**"GREEN"** means attaching well



## HMP-UNI350

### FEATURES:

- No drilling or screwing necessary
- Quickly installed on virtually or horizontally with no tools Provides assistance to the home or when traveling.
- Ideal for children while showering
- Support assist to over 150lbs
- Suction power Detective Device to prevent any unexpected injuries happened
- Large release levers make installing and removing easily.
- Anti - slide TPR over mold grip

### SPECIFICATIONS:

- Overall Length : 350 mm
- Suction Cup Diameter : 93 mm
- Constructed of heavy duty ABS plastic.
- Color: White or light blue

### Suction Power Detective Device Multinational Patent available

- Timely detect any loosen Suction Cup mounted on the surface in bathroom
- To prevent users against accidents caused
- Each of units includes a probe indicating member visible via the window
- Indicating window shows different colors

BATH & TOILET AIDS



### SPECIFICATIONS:

- Overall Length : 500 mm
- Suction Cup Diameter : 93 mm
- Color: White or light blue



**"RED"** is a warning signal request re-mount again



**"GREEN"** means attaching well

## HMP-UNI500

### FEATURES:


- No drilling or screwing necessary
- Quickly installed on virtually or horizontally with no tools Provides assistance to the home or when traveling.
- Ideal for children while showering
- Support assist to over 150lbs
- Suction power Detective Device to prevent any unexpected injuries happened
- Large release levers make installing and removing easily.
- Additional anti-slide TPR pads under the grip

## FEATURES:

1. No drilling or screwing necessary
2. Quickly installed on virtually or horizontally with no tools
3. Provides assistance to the home or when traveling.
4. Ideal for children while showering
5. Support assist to over 150lbs
6. Suction Power Detective device to prevent any unexpected injuries happened.
7. Large release levers make installing and removing easily.
8. Parent-Taiwan, Mainland, Japan, US, Germany, Spain..
9. CE marking.

**HMP-UNI400X**


**HMP-UNI500X**



# Suction Bathtub Handle

**HMP-UNI400S**

**HMP-UNI500S**



## SPECIFICATIONS:

- Overall Length
  1. HMP-UNI500S: 547 mm
  2. HMP-UNI400S: 445 mm
  3. HMP-UNI500X: 511 mm
  4. HMP-UNI400X: 409 mm
- Suction Cup Diameter: 93 mm
- Constructed of heavy duty ABS plastic.
- Color: As per requirement

Multi-national Patent  
Suction Indicator Device

- **“GREEN”** means attaching well
- **“RED”** is a warning signal request re-mount again

## Suction Bathtub Handle

### FEATURES:

- No drilling or screwing necessary
- Quickly and easily assembly with no tools
- Any type of composition/connection according to the size of surface in bathroom
- Innovative steering head suction cup unit allows pipe in different direction installed.
- Suction power Detective Device to prevent any unexpected injuries happened
- Side slide lever with locker to keep suction cup from removing easily.
- Anti-slide TPR pads under the grip
- Support assist to over 150lbs
- Other accessories: Shower holder, tower hanger, soap container etc. could be fitted up with UNI Clip

### SPECIFICATIONS:

- Overall Length:  
HMP-UNI600D: 500 mm  
HMP-UNI800L : 890 mm
- Suction Cup Diameter: 93 mm
- Color: White base and light blue C-clip and side lever

BATH & TOILET AIDS


### Suction Power Detective Device Multinational Patent available

- Timely detect any loosen Suction Cup mounted on the surface in bathroom
- To prevent users against accidents caused
- Each of units includes a probe indicating member visible via the window
- Indicating window shows different colors
- **"GREEN"** means attaching well
- **"RED"** is a warning signal request re-mount again

## HMP-UNI600D



## HMP-UNI800L



Innovative steering head



allows pipe connection in different direction installed.



Position fixed with screw nut



## FEATURE

- \* Multi-functional shower chair designed to combine the tub bar, shower chair and swivel chair as one with Max. Weight of 300lbs.
- \* Tool-less assembly on the top of wall of bathtub.
- \* Safety handle provides support and balance when entering and exiting bathtub.
- \* Swivel seat provides easily getting in and out of bathtub for the elder or disable people.

## HMP-SC301-A Swivel Shower Chair /Long Tub Bar




### Specifications:

Product size: 500 (L) \* 312 (W) \* 500 (H) mm

### Materials:

- 1.Seat: PE
- 2.Matress:TPE
- 3.Handle:TPR
- 4.Others: A3 steel

## HMP-SC301-B Swivel Shower Chair / Short Tub Bar



### Specifications:

Product size: 535 (L) \* 312 (W) \* 360 (H) mm

### Materials:

- 1.Seat: PE
- 2.Matress:TPE
- 3.Handle:TPR
- 1.Others: A3 steel

## HMP-SC304-A Swivel Seat



### Product size:

L: 505 mm  
W: 346 mm  
H: 405 ~ 532 mm

### Materials:

1. Matress:TPE
2. Handle:TPR
3. Leg: Aluminum

## HMP-SC304-B Swivel Seat (Short Tub Bar)



### Product size:

L: 505 mm  
W: 346 mm  
H: 513~ 638 mm

### Materials:

1. Matress:TPE
2. Handle:TPR
3. Leg: Aluminum

## HMP-SC302 5 Increment Height Adjustment Tub Bar



### Product size:

L: 218 mm  
W:170 mm  
H: 568~ 483 mm

### Materials:

1. Matress:TPE
2. Handle:TPR
3. Others: A3 steel

## HMP-SC303 Height Fixed Tub Bar



### Product size:

L: 203 mm  
W: 170 mm  
H: 524 mm

### Materials:

1. Matress:TPE
2. Handle:TPR
3. Others: A3 steel



- FEATURE**
- \* No drilling or screwing necessary, quickly installed in the bath tub.
  - \* Most innovative way to benefit fitness and health by exercise or rehab 10 minutes while bathing everyday
  - \* Help to the blood circulation/ beauty treatment/ body shape/ slim and muscle exercise.
  - \* Prevent sports injury which occurred easily on land exercise due to the water floating power and less-friction in the bath
  - \* Ideal for the aged people's rehabilitation or exercise therapy
  - \* \* Patented safety controller could timely detect the rail's adhesiveness:
- “**GREEN**” means the hand rail has been correctly mounted.
- “**RED**” means the hand rail is loose and needs to be reapplied.

- SPECIFICATIONS:**
1. Overall Length: 430mm
  2. Big S-Cup Diameter: 93mm
  3. Small S-cup Diameter: 50mm

- MATERIALS:**
1. ABS hand rail
  2. PVC Big Suction Cup
  3. TPR Small Suction Cup

