

## TRANSFER WOOD BOARD

**HMP-TS30910**

Easy for user to transfer between bed and wheelchair or...

Material: Wood  
Size: L =610~730 mm  
W=260 mm  
H =15 mm

Package:  
5PCS/CTN  
Dimension:740\*360\*90mm  
N.W.: 7.7 KGS  
G.W.: 8.8 KGS



Front



Back

## WOOD FOOT STOOL

**HMP-TS30920**

Material: Wood  
Size: L =392 mm  
W=302 mm  
H =105 mm

Package:  
6PCS/CTN  
Dimension:615\*405\*325mm  
N.W.: 15 KGS  
G.W.: 17.4 KGS



Front



Back

## FOLDING SHOWER CHAIR

**HMP-5000F**

Size:SEAT HEIGHT=460~505 mm  
SEAT WIDTH=405 mm  
SEAT DEPTH =305 mm

Package:  
4PCS/CTN  
Dimension:525\*435\*415mm  
N.W.: 9.6 KGS  
G.W.: 11 KGS



# BED STICK WITH WOOD/PLASTIC BOARD

Material: Wood and Steel

Wood Size:

L = 160 mm

W = 255~325 mm

Wood Thickness = 15 mm

Rail Size:

H 420-510mm; W 230mm

Package:

6PCS/CTN

Dimension: 640\*370\*400mm

N.W.: 13 KGS

G.W.: 15.5 KGS

## HMP-BS100



## HMP-BS200



Home Medical Products Inc.

[Www.homemedical.com.tw](http://www.homemedical.com.tw)