



Home Medical Products Inc.

FOLDING WALKER



Economic Folding Walker



HMP-6001

HMP-6001AW5

Folding Walker (Junior)



HMP-6001J

HMP-6001AW5

Reciprocal Motion Folding Walker



HMP-6001R

WALKER



HMP-6002

Deluxe Folding Walker


Bearing Capacity: 250 lbs
 Slightly flared legs enhance stability.
 Single push button ,easy operation.
 Adjustable height : 1" increments
 from 31" - 36" (ground to handgrip)
 Surface Processing: Shiny or Satin
 finish available
 Color: Silver, Bronze and Black available.



HMP-6002R

Reciprocal walker

Bearing Capacity: 250 lbs
 Slightly flared legs enhance stability.
 Single push button ,easy operation.
 Adjustable height : 1" increments
 from 31" - 36" (ground to handgrip)
 Surface Processing: Shiny or Satin
 finish available
 Color: Silver, Bronze and Black available.



HMP-6003

Reciprocal Motion Folding Walker

Bearing Capacity: 250 lbs
 Slightly flared legs enhance stability.
 Single push button ,easy operation.
 Adjustable height : 1" increments
 from 31" - 36" (ground to handgrip)
 Surface Processing: Shiny or Satin
 finish available
 Color: Silver, Bronze and Black available.



HMP-6004

Economic Fixed Walker

Bearing Capacity: 300 lbs.
 Slightly flared legs enhance stability.
 Adjustable height : 1" increments
 from 30" - 34" (ground to handgrip)
 Extra light weight: 2 kg
 Wheels able to be added on front legs.



Home Medical Products Inc.

WALKER



HMP-6005

Folding Side Walker

Bearing Capacity: 250 lbs.
The front legs slightly bended for touching ground vertically.
Adjustable height : 1" increments from 29" -33" (ground to handgrip)
Wheels able to be added on front legs.



HMP-6006A

Dual Release Button Folding Walker (Adult)

Bearing Capacity: 300 lbs
Slightly flared legs enhance stability.
Single push button ,easy operation.
Adjustable height : 1" increments from 32" - 39" (ground to handgrip)
Surface Processing: Shiny or Satin finish available
Color: Silver, Bronze and Black available.



HMP-6006MB

folding walker w/5" swivel castors

Bearing Capacity: 300 lbs
Slightly flared legs enhance stability.
Single push button ,easy operation.
Adjustable height : 1" increments from 34" - 38" (ground to handgrip)
Surface Processing: Shiny or Satin finish available
Color: Silver, Bronze and Black available.



HMP-6006X

Dual Release Button Folding Walker (Extra Wide)

Bearing Capacity: 400 lbs.
Slightly flared legs enhance stability.
Adjustable height : 1" increments from 30" - 34" (ground to handgrip)
Extra Wide: Width 22" between hand grips
Wheels able to be added on front legs.



Home Medical Products Inc.

WALKER



HMP-6008

Walker with Double Decker Handles

Bearing Capacity: 300 lbs.
Adjustable height : 1" increments from 32" -35" (ground to handgrip)
Surface Processing: Shiny or Satin finish available
Color: Silver, Bronze and Black available.



HMP-6012

Walker with single H bar

1" aluminum armrest tube with silver finish.
Grey single H shape bar with powder coating finish.
Adjustable height : 1" increments from 28" - 35" (ground to handgrip)
Surface Processing: Shiny or Satin finish available
Color: Silver, Bronze and Black available.



HMP-6013

Dual Purposes Folding Walker

Bearing Capacity: 300 lbs
Slightly flared legs enhance stability.
Single push button ,easy operation.
Adjustable height : 1" increments from 28" - 35" (ground to handgrip)
Surface Processing: Shiny or Satin finish available
Color: Silver, Bronze and Black available.
Dual purposes: Users are free to choose to move on straight way or reciprocally.

